


How to make a tech QA note.


Tech QA notes can be created in the patient section of Clario. This is located in the bottom left of the following image, located where communication notes are created.

The screenshot displays the Clario patient interface. At the top, it shows patient information: Last Name: *Test* (241 RVU:37.95 WU:239.00), TEST, MARY; Patient: TEST, SMOKER SUZIE; Site: CHW; and Ordering: ANDREW SPR. Below this, patient details are listed: MRN: 005590357CHWW, Gender: Female, Age: 73y, DOB: 01/03/1946, MPI: 217364, SSN, Home #, Work #, and Mobile #. The main section contains a report for a CT CHEST FOR LUNG CANCER SCREENING. The report includes sections for EXAM, INDICATION, TECHNIQUE, COMPARISON, FINDINGS, and IMPRESSION. The FINDINGS section lists: LOWER NECK AND AXILLA: Unremarkable; HEART AND VESSELS: Unremarkable; MEDIASTINUM AND HILA: Unremarkable; LUNG AND AIRWAYS: There are multiple bilateral small pulmonary nodules which have a tree-in-bud morphology; PLEURA: Unremarkable; BONES AND CHEST WALL: Unremarkable; UPPER ABDOMEN: Unremarkable. The IMPRESSION section states: There are bilateral tree-in-bud infiltrates which can be seen with bronchiolitis. Correlation with patient's current respiratory symptoms recommended. These numerous nodules have the potential of obscuring true pulmonary nodules. At the bottom left, a 'Notes' section is visible with a toolbar containing icons for communication and a 'QA Note Icon' highlighted by a red arrow. The QA Note Icon is a small circular icon with a white background and a black outline, located in the bottom left corner of the interface.

Select the Icon that looks like a First Place ribbon.




Once the QA Note pop-up opens, modality will be auto entered.
A section for your note is at the bottom.



3

Select from one of the three choices of either documentation error, imaging protocol error, or positioning error.




Technologist QA:

Rating:

Action Taken:

<Note Message>

Once your QA note is submitted you will see your note in the same section below the create note icons. This note will be exported in a report for next actions.



Notes

TECH QA NOTE - Positioning error Pending

Error in positioning

Carr, Matt, 02/21/2019 9:27 AM

If you have any questions about how to create a QA note, please call your CSA team at.

1-800-764-4563

