

MR Enterography

Fluids only 12 hours prior

NPO 4 hours prior

FOV: Start at pubic symphysis....48cm superior...if miss a little liver, that's ok

- Cor Cine
 - 15 images / table position
- SSFSE 40mm thick overlap 10mm
- Give 1mg Glucagon IM
- Ax FIESTA 5mm gap 1 265x224
- Cor FIESTA 3mm gap 1 256x224
- As SSFSE 4/1 256x224
- Ax SSFSE Fat Sat 4/1
- Cor SSFSE FAT3/1
- Lava Cor and Ax 2.5/1 256x224
 - 70 sec delay
 - 2cc/sec

Full protocol to follow